

LABORATORY BLENDERS

1 & 1.2 Liter Blenders

COMMON FEATURES:

Two-speed, push-button operation

Base: die-cast, epoxy-coated motor housing

Mechanical timer

.5 HP motor



SKU	Technical Information	Approvals	Plug
7010G	120V, 60 Hz	cTUVus	B ☹️
8010G	230V, 50 Hz	RoHS	F ☹️
8010EG	230V, 50 Hz	CE, RoHS	F ☹️

- 1.2-liter, heat-resistant glass container with handle and two-piece vinyl and styrene lid

SKU	Technical Information	Approvals	Plug
7010S	120V, 60 Hz	cTUVus	B ☹️
8010S	230V, 50 Hz	RoHS	F ☹️
8010ES	230V, 50 Hz	CE, RoHS	F ☹️

- 1-liter, stainless steel container with handle and two-piece vinyl and styrene lid